



Fluro Image of Actual Epidural Placement
courtesy of Claudio Feler, M.D.

14g RX™ Coudé®

The patented RX™ Coudé® Epidural Needle's rear heel is rounded to resist catheter/lead shearing during placement.



Product Options

Lengths: 4" Catalog # 105-1314

6" Catalog # 105-1614

Features

- ▶ Coudé® (Curved)
- ▶ Radiopaque
- ▶ Depth Markings & Attachable Wing
- ▶ Printed arrows on hub indicate direction of curve
- ▶ Easily Removable Stylet



New! RX²™

14g RX²™ Blunt Stylet



Extended Blunt Stylet

www.epimedpain.com

Catalog # 107-1314

The RX²™ blunt obturating stylet extends from the tip to push tissue structures away from the cutting edge during needle rotation. Simply place the locking obturating stylet into the RX™ Coudé® Needle after epidural placement is confirmed.